

AirFit™ for Her mask series

With ResMed's new AirFit for Her mask series, it's easier than ever to find the right fit for you. Choose from the AirFit P10 for Her nasal pillows mask, the AirFit N10 for Her nasal mask or the AirFit F10 for Her full face mask for a lightweight, stylish solution designed to make therapy comfortable and easy – no matter which one you choose.



sleepvantage
your sleep community

Membership. Support. Community.

Hundreds of thousands of people, like you, require treatment for their sleep apnea. It's time to get connected through our sleepvantage member support program.

Join sleepvantage to receive the support you need and enjoy a range of member benefits. It's absolutely free!

Join today

🌐 sleepvantage.com.au

☎ 1300 305 705 (Aus) or 0800 737 633 (NZ)

✉ info@sleepvantage.com.au



ResMed Air Solutions

great sleep



a new beginning
in great sleep



- 1 Valipour A et al. Gender-related differences in symptoms of patients with suspected breathing disorders in sleep: a clinical population study using the sleep disorders questionnaire. *Sleep* 2007; 30(3): 312-9
- 2 Callop N, Cassel DK. Snoring and sleep disordered breathing. In Lee-Chiong T Jr, Sateia M, Carskadon M. Editor. *Sleep Medicine*. Philadelphia: Hanley & Belfus, 2002: 349-355
- 3 Mihai V, Rusu G, Mihaescu T. Demographic clinical and polysomnographic differences between men and women. *Pneumologia* 2010; 59(2): 64-7
- 4 Lin CM, Davidson TM, Ancoli-Israel S. Gender differences in obstructive sleep apnea and treatment implications. *Sleep Med Rev* 2008; 12(6):481-96

ResMed Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

See www.ResMed.com for ResMed locations worldwide. AirFit, AirSense, AutoRamp, AutoSet, ClimateLine and HumidAir are trademarks of ResMed Ltd. AirFit, AutoSet and ClimateLine are registered in the U.S. Patent and Trademark Office. Specifications may change without notice. Product availability may vary across regions. ©2014 ResMed Ltd. 1017611/1 2014-07

ResMed.com/AirSolutions

ResMed.com/AirSolutions

Research shows that sleep apnea in women is different to sleep apnea in men.¹ So it makes sense that your treatment should be different too. ResMed Air Solutions introduces the first complete sleep therapy solution tailored for women — combining cutting-edge technology with stylish masks and personalised support to give you everything you need for a great night's sleep.



Unique therapy for women

Studies show that women with obstructive sleep apnea are less likely to snore, and have more flow limitations and frequent arousals during sleep than men.^{2,3,4} To address these important differences, ResMed has created the sensitive AutoSet™ for Her algorithm which detects flow limitation as early as possible, and responds by applying pressure changes steadily and comfortably throughout the night.

Out-of-the-box comfort

The AirSense™ 10 AutoSet for Her has been designed to make therapy easy. With no complicated settings or menus to navigate, simply press Start and you're ready to go. Each device comes with the new HumidAir™ heated humidifier built-in, so you can enjoy the benefits of humidification in one easy-to-use system. And when used with the ClimateLineAir™ heated tube, your device is designed to automatically deliver the optimal temperature and humidity, for a more comfortable experience.

A comfortable start

The new AutoRamp™ feature with sleep onset detection helps make therapy comfortable from the moment you turn on your device. AutoRamp works by delivering a low pressure to help you fall asleep with ease. Once it detects you're asleep, it comfortably increases the pressure to ensure your prescribed level is delivered the moment you need it.

A restful night

Every device has a built-in ambient light sensor that adjusts the brightness of the screen based on the amount of light in the room. So if you get up during the night and need to pause your therapy, the screen will stay dim so it doesn't disrupt your bed partner. The screen's back light also turns off automatically to help you fall back to sleep without disturbance. And with ResMed's Easy-Breathe motor built into each device, you and your bed partner can enjoy whisper-quiet sleep, every night.



Find out more about
ResMed Air Solutions at
[ResMed.com/AirSolutions](https://www.resmed.com/AirSolutions)