

## End of support life – Stellar Type 1

### Details

Product support for Stellar Type 1 devices has ended, or will soon be ending. The end of support dates are as follows:

Product code	Product	End of support date
24140	Stellar 150 Americas	14 <sup>th</sup> June 2018
24147	Stellar 150 APAC	20 <sup>th</sup> June 2019
24156	Stellar 100 APAC	30 <sup>th</sup> May 2018
24157	Stellar 100 Americas	4 <sup>th</sup> February 2018

The following product codes had both Type 1 and Type 2 variants. The Type 1 devices prior to the serial number break listed are no longer supported.

Product code	Product	Serial number break	End of support date
24141	Stellar 150 Germany	20122338703	22 <sup>nd</sup> August 2017
24142	Stellar 150 France	20122338763	21 <sup>st</sup> August 2017
24143	Stellar 150 EUR G1	20122329483	22 <sup>nd</sup> August 2017
24144	Stellar 150 EUR G2	20122338885	8 <sup>th</sup> August 2017
24145	Stellar 150 UK	20122338887	30 <sup>th</sup> July 2017
24151	Stellar 150 Germany	20122338765	7 <sup>th</sup> August 2017
24152	Stellar 150 France	20122338825	2 <sup>nd</sup> August 2017
24153	Stellar 150 EUR G1	20122338889	18 <sup>th</sup> May 2017
24154	Stellar 150 EUR G2	20122338891	25 <sup>th</sup> December 2016
24155	Stellar 150 UK	20122338893	24 <sup>th</sup> August 2017

While Stellar 1.0 devices may still function as intended, the availability of spare parts and repair service for these devices cannot be guaranteed beyond these dates.

Type 1 units can be distinguished from Type 2 units as they do not have Remote Alarm support.



Type 1 unit



Type 2 unit with support for Remote Alarm

**Note:** Support remains for Type 2 Stellar devices.