



Data management/Device compatibility list

English

This supplement specifies the devices compatible with ResMed data management products. For full details regarding the correct use of these devices, see the relevant User Guides.

Note: Not all products are available in all regions.

	ResScan™ 6.1	EasyCare Tx / Tx Link 9.0	AirView™	RCM	RCMH*
AirMini™			✓		
AirSense 10 / AirCurve 10 Series					
AirSense™ 10 AutoSet / AirSense™ 10 AutoSet for Her	✓	✓	✓		
AirSense™ 10 Elite	✓	✓	✓		
AirCurve™ 10 CS PaceWave™	✓	✓	✓		
AirCurve™ 10 CS-A PaceWave™	✓	✓	✓		
AirCurve™ 10 S	✓	✓	✓		
AirCurve™ 10 ST	✓	✓	✓		
AirCurve™ 10 VAuto	✓	✓	✓		
Lumis Series					
Lumis™ 100 VPAP S	✓	✓	✓		
Lumis™ 100 VPAP ST	✓	✓	✓		
Lumis™ 150 VPAP ST	✓	✓	✓		
Lumis™ 100 VPAP ST-A	✓	✓	✓		
Lumis™ 150 VPAP ST-A	✓	✓	✓		
Lumis™ ST EU 30cmH2O PAC			✓		
S9 Series					
Auto 25 (S9)	✓	✓			
AutoSet CS-A (S9 PaceWave™)	✓	✓			
AutoSet CS™ (S9 PaceWave™)	✓	✓	✓		
AutoSet CS™ (S9)	✓	✓	✓		
S9 AutoSet™	✓	✓	✓		
S9 Elite™	✓	✓	✓		
S9 Escape™	✓	✓			
S9 Escape™ Auto	✓	✓			
S9 VPAP™ Tx		✓			
V-Auto (S9)	✓	✓	✓		
VPAP™ Adapt (S9 PaceWave™)	✓	✓			
VPAP™ Adapt (S9)	✓	✓			
VPAP™ ST (S9 iVAPS)	✓	✓	✓		
VPAP™ ST-A (S9 iVAPS)	✓	✓			
VPAP™ ST-A (S9)	✓	✓	✓		
VPAP™ S (S9)	✓	✓			
VPAP™ ST (S9)	✓	✓	✓		
S8 Series					
S8 AutoScore™	✓	✓			
S8 AutoScore™ II	✓	✓			
S8 AutoSet Spirit™	✓	✓			
S8 AutoSet Spirit™ II	✓	✓			
S8 AutoSet Vantage™	✓	✓			



	ResScan™ 6.1	EasyCare Tx / Tx Link 9.0	AirView™	RCM	RCMH*
S8 AutoSet™ II	✓	✓			
S8 Elite™	✓				
S8 Elite™ II	✓	✓			
S8™ Auto 25	✓	✓			
VPAP™ IV Series	✓	✓			
S7 Series					
AutoSet CS™ 2	✓	✓			
S7 AutoSet Spirit™	✓				
S7 Elite™	✓				
VPAP™ Adapt SV	✓				
VPAP™ II Series	✓				
VPAP™ III Series	✓	✓			
VPAP™ III ST-A with QuickNav	✓	✓			
Ventilator					
Astral™ 100	✓		✓	✓	✓
Astral™ 150	✓		✓	✓	✓
Stellar™ 100	✓	✓	✓	✓	
Stellar™ 130	✓	✓	✓	✓	
Stellar™ 150	✓	✓	✓	✓	
Elisée™ 150	✓				
VS III™	✓				
Home sleep testing device					
ApneaLink™ Air			✓		

RCM – ResMed Connectivity Module

RCMH – ResMed Connectivity Module Hospital

*For a list of available data for use by the EMR system, please contact your ResMed representative.

 **MANUFACTURER** ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirMini, AirSense, AirCurve, AirView, Astral, ApneaLink, Lumis, S8, S9, Stellar, VPAP, AutoSet, S9 Escape, S9 Elite, S9 Escape Auto, S8 Elite, S8 AutoSet Vantage, PaceWave, ResScan, VPAP Tx, S8 AutoSet, S9 AutoSet, S9 Auto 25, VS III, Elisée, VPAP Adapt SV, AutoSet CS, S8 AutoScore, S8 AutoSet Spirit and VPAP are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip.
1018030/19 2020-08